

## KERi Complete Nursing Skills Manikin, light skin

### Product Information

Weight	13.5 kg
Brand	Nasco
Item No.	W44075 [1013742]

- This product carries a 5 year warranty



### Description

KERi™ has a non-age-specific appearance. Complete lifelike range of motion, realistic patient positioning, and non-pinching joints - moves like a person.

- Lightweight - approximately 28 lbs.
- Full-size adult manikin - measures 58"
- Overall female appearance
- Converts to male with removal of wig and attachment of male genitals.

### The following nursing and patient care procedures can be performed:

- Prostate examination (Stage B)
- Ostomy Care (Ileostomy and colostomy lavage and suctioning)
- Tracheostomy care (lavage and suctioning)
- Placement of Nasogastric tubes
- Oral and nasal lavage, gavage and suctioning
- Enema administration (female), pap smears and douching
- Intramuscular Injection (arms, thigh, buttock)
- Oral hygiene (removable upper and lower denture)
- Eye irrigation, bed baths, hair care (washing, combing)
- Patient positioning and transfer techniques, clothing changes
- Ear irrigation and hearing aid placement and removal
- Bandaging and wound dressing (including fingers and toes)
- Urinary Catheterization - Female and Male (Uncircumcised)

### Visual Inspection Features include:

- Decubitus Sacral Ulcer - Stage 1
- Dilated Pupil Comparison

- Cancerous Mole Comparison
- Reddened Skin Folds

**Complete Range of Motion include:**

- Trunk - Rotation, Hyperextension
- Neck - Rotation, Hyperextension, Lateral Flexion
- 

---

[https://www.3bscientific.com/no/keri-complete-nursing-skills-manikin-light-skin.p\\_153\\_4950.html](https://www.3bscientific.com/no/keri-complete-nursing-skills-manikin-light-skin.p_153_4950.html)

More product images

