

works.



WORK STOOLS AND WORK CHAIRS FOR SURGERY, CLINICS AND LABORATORIES.



GREINER

workstools.



MEDICAL WORK STOOLS AND CHAIRS – MADE TO MEASURE.

The variety of the GREINER work stools programme is as versatile as the working life in the medical field. Well known for their „Made in Germany“ experience, the stools have been supporting the work of medical professionals 24/7 around the world for years. The range of applications extends from examination and treatment up to the operating room. Both because of their strengths in anatomy, ergonomics, and functionality and because of accommodating GREINER medical chairs for specialist applications. In other words: day after day, GREINER helps to make the everyday work of healthcare better, easier, and healthier. Good seating in the best quality which facilitates processes while strengthening the back.



WE DESIGN QUALITY.

With innovative features, GREINER has always managed to show that there is a difference in work stools. This includes the fulfilment of hygienic and technical requirements for treatment as well as the varied selection of heights for sitting and upholstery, height adjustment with a 360° manual ring or foot release of the gas spring as well as the choice of base part and castors. Work stools and chairs of high-class quality that retain a long working life are produced at the in-house facilities including a foaming and sewing department. The cover in premium material is extremely comfortable, which noticeably improves the sitting comfort. Of course, all work stools meet the clinical requirements for stability, safety and hygiene, while sustainability in production and permanent spare parts supply are guaranteed.

Medical work stools and work chairs.

- ◆ A wide selection of models with hand or foot release, fixed or mobile version, with swivelling and height-adjustable backrest
- ◆ Great comfort and fulfilment of requirements for sterility through innovative height adjustment: ergonomic ring release by hand or a ring-shaped base plate that can be released by foot
- ◆ Three different sitting heights to fit body size, seat cushions with a diameter of 36 cm or 40 cm and an angled anatomically-shaped seat edge included
- ◆ Models with a full-foam seat and backrest for special hygiene requirements
- ◆ Safety castors keep the work stool in position when it is not being moved
- ◆ Choice of cover materials compact and premium: see colour pattern card



A DESIGN TO VITALISE YOUR POSITION.

Sit like a pro so you can work like a pro. Anatomy, ergonomics and functionality characterise the product design of GREINER's work stools and chairs. Modern sitting comfort not only stands for comfortable upholstery and backrests or smooth-running castors – every functional detail also contributes to an active and ergonomic sitting position.



Versatile work stool that suits every treatment chair




Large backrest front face



Large backrest from the rear



Foot release with foot ring



Stainless steel round plate



Fatigue-free operation thanks to support for arms or torso



Surgery-compliant equipment EL 11 with electrical conductivity



"Black" conductive cover, EL 11



Electrically conductive safety castors, EL 11



Crescent-shaped backrest
Combined back/armrest, swivelling and height adjustable

Anatomical sitting cushion
Perfect ergonomics for examination and treatment

Manual release
Height adjustment with chrome-plated handle for smooth and precise movement

Tailored sitting height
Column options with different heights can be selected

Perfect matching base parts
Chromed foot ring here, see below for other options

Safety castors
Smooth-running castors with stop function when the load changes

Base parts according to choice
5-star base made of polished aluminum; for options see below

Doubtful if anyone will then say "All stools are the same." Let's take a close look at the GREINER work stools' design – from bottom to top: safety castors with/without central locking or sliders? Short, medium or long gas springs for different sitting positions? Height adjustment with manual ring release or "hands-free" by foot? Sitting cushion, padded and sewn by hand – round or anatomically shaped, cover with compact or premium fabric, or a full foam seat – super soft, with a hygienic, breathable surface? In addition, the backrest, whether small or large, swivelling as a back/armrest in the shape of a crescent moon, or completely without it? And finally the question of electrical conductivity of the cover, column and castors that is suitable for surgery? Any other questions? Time to have a look at the details.



Sitting cushion Ø 36 or 40 cm
Cover in compact or premium fabrics



360° manual ring release
Height adjustment all around from any sitting position

Chromed column
High-quality and durable surfaces

VARIETY IS KEY - THE BASE PARTS.

This is what free choice looks like: The base parts of the GREINER work stools and chairs are characterised by their glossy design, quality and seating comfort. Otherwise, they are rather versatile and tailored to your needs. The foot release developed by GREINER, for example, offers high clinical added value wherever sterile or contaminated gloves are used. Triggered by the round pedal plate attached above the 5-star base, it can be operated "hands-free" all around. All work stools have safety castors, on request as an electrically conductive version.



Manual ring release, 5-star base Manual ring release, foot ring Foot release with foot ring 5-star base with sliders Central arrest



NexGen design and manual ring release



Crescent-shaped full foam backrest



Modern design with high functionality



Breathable cushion

FULL FOAM – PUR(E) ADDED VALUE.

Being innovative means always re-evaluating products and their benefits against the backdrop of changing requirements or the development of advanced technologies and materials. And if a pioneering generation of materials proves to be ready for the market, they can also be implemented in new products. With the new PUR integral (full) foam, a material from the group of polyurethanes for sitting cushions and backrests is now available for medical use, which in addition to a soft and breathable core has an increasingly firmer texture towards the outside.

Full foam seats guarantee a high level of sitting comfort, a pleasant feeling on the surface, are robust in handling and free of seams – making them very easy to clean and disinfect. The material is light and facilitates "dynamic" sitting where you can actively change your position. Easy-care and washable, PUR is characterised by its durability against many types of cleaning agents and disinfectants.



Pure dynamic back

The crescent-shaped fixed backrest offers you a secure hold and can also be used as an armrest.

Pure innovation

Full foam guarantees a non-slip hold thanks to closed, shock and waterproof surfaces. In addition, it is characterised by its durability against many types of cleaning agents.

Pure sitting comfort

The seat with 85 mm thick PUR full foam is pleasantly soft, free of seams and feels particularly good when you have to sit on it for hours.



Pure ease

This work stool with foot release is made from special integral foam. It is pleasantly light. Its air ventilating texture guarantees a pleasant feeling when seated.





The examples shown here represent just a small selection out of the broad variety.

Let us offer you some inspiration: Here is a small selection from the GREINER work stools and chairs' product range. Stools without backrest, chairs with backrest in different designs. Height adjustment by hand or foot. Base parts and sitting heights of your choice. Upholstery materials in a wide variety of colours with compact or premium fabric. Or EL 11 equipment that is suitable for surgeries with electrically conductive components. The choice is yours.



WORK STOOL 54 - 73 cm.
Item code 3008305
5-star base polished, sitting cushion Ø 36 cm, cover in compact „Petrol“ C435, 360° manual ring release, safety castors.



WORK STOOL 54 - 73 cm.
Item code 4108305 RG 01 PREMIUM HO
5-star base polished with sliders, sitting cushion Ø 40 cm, cover in premium „Aqua“ P700, 360° manual ring release.



WORK STOOL 54 - 73 cm.
Item code 3008305
5-star base polished, sitting cushion Ø 36 cm, cover in compact „Mint“ C450, 360° manual ring release, safety castors.



WORK STOOL 51 - 71 cm.
Item code 4105305 EL 11
Foot release, sitting cushion Ø 40 cm, cover „Black“ EL 11, gas spring column and safety castors overall electrically conductive.



WORK STOOL 54 - 73 cm.
Item code 3009105
Chromed foot ring, sitting cushion Ø 36 cm, cover in compact „Citron“ C500, 360° manual ring release, safety castors.



WORK STOOL 40 - 52 cm.
Item code 4108105 PREMIUM HO
5-star base polished, sitting cushion Ø 40 cm, cover in premium „Aqua“ P700, 360° manual ring release, safety castors.



WORK STOOL 54 - 73 cm.
Item code 4129105
Chromed foot ring, anatomically shaped seat edge Ø 40 cm, cover in compact „Cyan“ C422, 360° manual ring release, safety castors.



WORK STOOL 46 - 58 cm.
Item code 3008205 CENTRAL
5-star base polished with central arrest, sitting cushion Ø 36 cm, cover in compact „Orange“ C530, 360° manual ring release.



WORK CHAIR 54 - 73 cm.
Item code 4019105
Chromed foot ring, full foam sitting cushion and crescent-shaped backrest, „Black“, 360° manual ring release, safety castors.



WORK CHAIR 40 - 52 cm.
Item code 3018105
5-star base polished, small backrest height adjustable, sitting cushion Ø 36 cm, cover in compact „Stone“ C214, 360° manual ring release, safety castors.



WORK CHAIR 51 - 71 cm.
Item code 4015305
Foot release, full foam sitting cushion and crescent-shaped backrest, „Black“, safety castors.



WORK CHAIR 51 - 71 cm.
Item code 3045305
Foot release, backrest height adjustable, sitting cushion Ø 40 cm, cover in compact „Atoll“ C706, safety castors.



WORK CHAIR 54 - 73 cm.
Item code 4185105 PREMIUM AS
Chromed foot ring, swivelling crescent-shaped backrest height adjustable, anatomically shaped seat edge Ø 40 cm, cover in premium „Navy“ + „Aqua“ P710/700, manual height adjustment, safety castors.



WORK CHAIR 51 - 71 cm.
Item code 4185305 EL 11
Foot release, swivelling crescent-shaped backrest height adjustable, anatomically shaped seat edge Ø 40 cm, cover „Black“ EL 11, gas spring column and safety castors overall electrically conductive.




OFFICE CHAIR 54 - 73 cm.
Item code 3054600 PREMIUM BS
5-star base polished, backrest adjustable in height and depth with dumping function, without armrest, cover in premium „Navy“ P710, safety castors.



OFFICE CHAIR 54 - 73 cm.
Item code 3054600
5-star base polished, backrest adjustable in height and depth with dumping function, without armrest, cover in compact „Mint“ C450, safety castors.

WORK STOOLS AND WORK CHAIRS FOR ERGONOMIC SEATING.

Ergonomic sitting comfort that will convince you: the wide range of work stools and chairs with or without backrest boosts your performance when examining and treating patients. GREINER offers you a product range that is processed at the highest quality, which you can select and adapt perfectly with a variety of practical features according to your requirements. In terms of shape and colour, you will definitely find the option that will best fit your environment. In addition, all upholstery fabrics are biocompatible, heat-resistant, disinfectant-resistant, urine and blood-resistant, scratch-resistant, dirt-resistant and durable. The fabric EL 11 „Black“ that is suitable for surgeries is made of carbon fibre material and is therefore electrically conductive.



Colours and fabrics
GREINER'S cover material is available in two options: "compact", which is characterised by comfort and plenty of options in 20 colours, as well as "premium", the innovative line that sets new standards with its super soft and, at the same time, far more durable surface in 12 colours.

AT A GLANCE.

An overview of Greiner work stools and chairs: Choose what you like from sitting cushions, backrests and a large range of sitting heights. The hand or foot release is also a criterion, as is the foot ring. All work stools have a highly polished aluminium 5-star base, a chrome-plated column and safety castors that keep the chair in position at standstill. Electrically conductive versions are also available for almost all models.



WORK STOOLS

WORK STOOL WITH MANUAL RING RELEASE, UPHOLSTERY Ø 36 CM

	Item code
• work stool 40-52 cm	3008105*
• work stool 46-58 cm	3008205*
• work stool 54-73 cm	3008305*
• work stool 46-58 cm with foot ring	3009005
• work stool 54-73 cm with foot ring	3009105

WORK CHAIR WITH SMALL BACKREST HEIGHT ADJUSTABLE, MANUAL RING RELEASE, UPHOLSTERY Ø 36 CM

• work chair 40-52 cm	3018105*
• work chair 46-58 cm	3018205*
• work chair 54-73 cm	3018305*
• work chair 46-58 cm with foot ring	3019005
• work chair 54-73 cm with foot ring	3019105

WORK STOOL WITH MANUAL RING RELEASE, FULL FOAM SEAT Ø 36 CM

• work stool 40-52 cm	4008105*
• work stool 46-58 cm	4008205*
• work stool 54-73 cm	4008305*
• work stool 46-58 cm with foot ring	4009005
• work stool 54-73 cm with foot ring	4009105

WORK STOOL WITH MANUAL RING RELEASE, UPHOLSTERY Ø 40 CM

• work stool 40-52 cm	4108105*
• work stool 46-58 cm	4108205*
• work stool 54-73 cm	4108305*
• work stool 46-58 cm with foot ring	4109005
• work stool 54-73 cm with foot ring	4109105

WORK CHAIR WITH BACKREST HEIGHT ADJUSTABLE, MANUAL RING RELEASE, UPHOLSTERY Ø 40 CM

• work chair 40-52 cm	3028105*
• work chair 46-58 cm	3028205*
• work chair 54-73 cm	3028305*
• work chair 46-58 cm with foot ring	3029005
• work chair 54-73 cm with foot ring	3029105

WORK STOOL WITH MANUAL RING RELEASE, ANATOMICALLY SHAPED SEAT Ø 40 CM

• work stool 40-52 cm	4128105*
• work stool 46-58 cm	4128205*
• work stool 54-73 cm	4128305*
• work stool 46-58 cm with foot ring	4129005
• work stool 54-73 cm with foot ring	4129105

WORK CHAIR WITH CRESCENT-SHAPED BACKREST ADJUSTABLE, MANUAL RING RELEASE, UPHOLSTERY Ø 40 CM

• work chair 40-52 cm	4154105*
• work chair 46-58 cm	4154205*
• work chair 54-73 cm	4154305*
• work chair 46-58 cm with foot ring	4155005
• work chair 54-73 cm with foot ring	4155105

WORK STOOLS

WORK CHAIR WITH CRESCENT-SHAPED BACKREST ADJUSTABLE, MANUAL RING RELEASE, ANATOMICALLY SHAPED SEAT Ø 40 CM

	Item code
• work chair 40-52 cm	4184105*
• work chair 46-58 cm	4184205*
• work chair 54-73 cm	4184305*
• work chair 46-58 cm with foot ring	4185005
• work chair 54-73 cm with foot ring	4185105

WORK CHAIR WITH MANUAL RING RELEASE, FULL FOAM SEAT Ø 36 CM AND CRESCENT-SHAPED BACKREST

• work chair 40-52 cm	4018105*
• work chair 46-58 cm	4018205*
• work chair 54-73 cm	4018305*
• work chair 46-58 cm with foot ring	4019005
• work chair 54-73 cm with foot ring	4019105

WORK STOOL WITH FOOT RELEASE, UPHOLSTERY Ø 40 CM

• work stool 43-56 cm	4105205
• work stool 51-71 cm	4105305

WORK STOOL WITH FOOT RELEASE, FULL FOAM SEAT Ø 36 CM

• work stool 43-56 cm	4005205
• work stool 51-71 cm	4005305

WORK CHAIR WITH BACKREST HEIGHT ADJUSTABLE, FOOT RELEASE, UPHOLSTERY Ø 40 CM

• work chair 43-56 cm	3045205
• work chair 51-71 cm	3045305

WORK CHAIR WITH CRESCENT-SHAPED BACKREST ADJUSTABLE, FOOT RELEASE, UPHOLSTERY Ø 40 CM

• work chair 43-56 cm	4155205
• work chair 51-71 cm	4155305

WORK STOOL WITH FOOT RELEASE, ANATOMICALLY SHAPED SEAT Ø 40 CM

• work stool 43-56 cm	4125205
• work stool 51-71 cm	4125305

WORK CHAIR WITH CRESCENT-SHAPED BACKREST ADJUSTABLE, FOOT RELEASE, ANATOMICALLY SHAPED SEAT Ø 40 CM

• work chair 43-56 cm	4185205
• work chair 51-71 cm	4185305

WORK CHAIR WITH FOOT RELEASE, FULL FOAM SEAT Ø 36 CM AND CRESCENT-SHAPED BACKREST

• work chair 43-56 cm	4015205
• work chair 51-71 cm	4015305

OFFICE CHAIR FOR CLINICS AND LABS, BACKREST HEIGHT ADJUSTABLE

• office chair without armrest 42-54cm	3054400*
• office chair without armrest 46-58cm	3054500*
• office chair without armrest 54-73cm	3054600*

* Central arrest available.

WORK STOOLS – accessories

• work stool in premium material	PREMIUM HO
• work chair in premium material	PREMIUM AS
• office chair in premium material	PREMIUM BS
• Electrically conductive covers (black) for work stools and chairs	EL 11
EL 11	<ul style="list-style-type: none"> – electrically conductive cover black – stainless steel round plate – plastic covered gas spring column black – electrically conductive safety castors
• stainless steel round plate	RONDE
• slider instead of castors	RG 01
• central arrest	CENTRAL*

Dimensions and technical data

• seat width & depth	Ø 36 – 40 cm
• sitting height	40 – 73 cm
• sitting height with foot release	43 – 71 cm
• full foam padding thickness	8,5 cm
• upholstered padding thickness	7,5 cm
• maximum load	150 kg
• weight	7 – 16 kg
• maximum width	67 cm
• castor diameter	50 mm



we design quality.

Agreeably health-promoting seating has motivated us since 1922. For decades we have been making high-quality, technically sophisticated couches within the medbest range. Our innovative devices for examination, dialysis, oncology, blood sampling, care and transport are all made in Germany. Our quality commitment is verified by an independently audited quality assurance system: GREINER is certified under DIN EN ISO 9001 and 13485.



GREINER GmbH
Wettestraße 1 · 74385 Pleidelsheim · Germany
T +49 7144 81120 · F +49 7144 8112 99
info@greiner-gmbh.de · www.greiner-gmbh.de